

Getting started with Scantraq

Thank you for purchasing Scantraq. If CD-ROM installation does not start automatically, run ..\install\setup32.exe. Administrator rights are required.

- **Select Scantraq and NOT FilmStar when installing for the first time.**
- **Install Scantraq before attaching the USB Security ID Module.**
- **There are seven versions of Scantraq and five supplemental servers; please select only versions compatible with your instrument(s).**

This table follows the same order as the installer.

1	Agilent (HP) 8453 {Mhp8453.exe}	<ul style="list-style-type: none"> • HPIB interface required
2	PE Lambda 2-950, IR 580-983, Servers {Mpe983.exe}	<ul style="list-style-type: none"> • If no URA, recommended for Lambda 650-950
2a	+ Ocean Optics Server	<ul style="list-style-type: none"> • OOIBase32 required; contact FTG Software for a web-link and installation password
2b	+ PE Spectrum 100 Server	<ul style="list-style-type: none"> • PE Spectrum Express required
2c	+ PE Spectrum GX/2000 Server	<ul style="list-style-type: none"> • PE Spectrum (16-bit) required
2d	+ PE Spectrum One Server	<ul style="list-style-type: none"> • PE Spectrum (32-bit) required
2e	+ PE Lambda XLS Server	<ul style="list-style-type: none"> • XLS Print Utility required
3	Varian Cary {Mcary5.exe}	<ul style="list-style-type: none"> • Minimum installation Cary WinUV (System + Scan application) required
4	Hitachi U-3210/3410 {Mh3410.exe}	<ul style="list-style-type: none"> • No support for newer instruments
5	PE Lambda 650-950 with URA {Mpe950.exe}	<ul style="list-style-type: none"> • Requires UV Winlab 6 with URA driver
6	PE Lambda 18/19 (AX PCI board required) {Mpe19AX.exe}	<ul style="list-style-type: none"> • Requires Lambda AX board and drivers; contact FTG Software • See accompanying document <i>Installing the PCI Card for Lambda 18/19 or PE Lambda 19 32-bit Installation</i> at ftgsoftware.com/issues.htm. • Original PE ISA board no longer supported
7	PE Lambda 1050 (URA/TDM support) {Mpe1050.exe}	<ul style="list-style-type: none"> • Requires UV Winlab 6 with URA and TDM drivers installed

After starting Scantraq, open a Scan Method file appropriate for the spectrophotometer you are using. When ready to actually acquire data, don't forget to turn off Test Mode in File..Configuration..Preferences.

Updates may be downloaded from our web site. Support e-mails include licensing and installation passwords (same passwords for CD-ROM and web download) at the bottom. **Bookmark <http://ftgsoftware.com/updates.htm> so you can easily get to our support pages.** Alternatively, if your Scantraq computer accesses the Internet, use the links in the Help menu. You can print additional manuals available in Adobe PDF format.

In addition to its Windows point and click interface, Scantraq includes development modules BASIC, Workbook, SQPlot, Database and Report Generator. The modules are described in on-line help and illustrated in supplied examples.

Contact FTG Software if you have any problems installing or running Scantraq. Also, please keep us informed if your e-mail address changes and if there other users who should receive support e-mails.

We suggest that you do not install Scantraq database support until you understand exactly what it does (see <http://ftgsoftware.com/database.htm>). Install FileMaker support (FM8Setup.exe) after installing both Scantraq and FileMaker Pro Version 8/9/10. If you don't have FileMaker Pro and want to see how it works (with some restricted functionality) install <http://ftgsoftware.com/files/fmrtsetup.exe>.